STI	UD	EN	IΤ				

EXPERIMENTS

EXPERIMENT 01 - HORIZONTAL PROJECTILE LAUNCH

Objectives:

- Recognizing the physical quantities involved in a horizontal projectile launch
- Checking the relationship between the physical quantities present in a horizontal launch.

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
<i>07</i>	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING
XX	62001226	01	Un.	DIGITAL TIMER AZB-30 USB (*)
XX	62001201	02	Un.	PHOTOELECTRIC SENSOR PGS-D10 (*)
XX	04002037	01	Un.	FLIGHT TIME SENSOR TFS-D10


(*) It does not accompany the product. It is sold separately.

Part I - Determining the horizontal launch velocity (v_0) using the reach measure (A) and the launch height (y_0)



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the launcher so that the projectile has space to move and fall on the floor, as shown in the figure.





- 2. Measure the launch height (y_0) in relation to the floor. The measurement shall be made from the lower part of the ball indicated at the exit of the cannon, to the floor, according to the figure.
- **3.** With a plumb bob, mark on a paper sheet (scotch taped on the floor) the position (origin of horizontal displacement). The plumb bob must match the vertical passing through the center of the ball.



- **4.** Place the plastic ball in the cannon and compress the spring to the second stage. Pull the trigger and observe where the ball touches the floor.
- **5.** At the dropping point of the ball place a paper sheet (scotch tape it as well) with a carbon paper over it.
- **6.** Repeat the launch 10 times and measure the reached distance range A using the measure tape. Complete the table 1 with the values.

Table 1
Reached distance (m)

1
2
3
4
5
Average value

7. What are the active forces in the sphere after the launch?



STUDENT	
---------	--

_								
8.	Classify the	projectile	motion	according	to the two	o directions	(vertical a	and horizontal).

9. Which equation must be used in each motion?


10. Use the horizontal launch equations and calculate the time in the air.

Part II - Determination of horizontal velocity of launch (v_o) using the measurement of projectile passage time by two sensors




Do not look directly at the cannon exit, as it may be charged!

1. Assemble the projectile launcher as shown in the figure (sensors are positioned approximately 5cm apart):



- **2.** Align correctly the sensors.
- **3.** Connect the sensors to the timer.
- **4.** Measure the distance between the two sensors centers.
- **5.** Turn on the timer and check if it has identified the sensors. Test the sensors by interrupting the infrared passage. An intermittent exclamation mark (!) should appear in the upper right corner of the screen indicating that the sensors are identified.
 - When you turn on the timer, it displays the STANDBY screen and should show:
 - STB
 - Vb: 8,5V Output Voltage.
 - Sns: 2- Number of sensors identified.
 - Bobbin: off indicating the coil is off.





6. Select function 1 by clicking FUNC. The display will show the function 1 F1. The display should read "[F1]" "Btn nS:2". In this function the timer will be using two sensors, and the chronometer will present two time measurements. (See information in the timer manual).



- **7.** Schedule the timer for data acquisition by starting the time counting on the first sensor and interrupting the second one.
- **8.** Press SETUP, press for a few seconds until the display shows CFG, use the key (↓) and adjust from Btn to Sns. Click START to end programming. On the nS: 1 screen, it indicates that the timer will only take one measure of time.



- **9.** As the objective is only to measure the time interval between the two sensors, collect the ball just after passing through the second sensor, using a cardboard box.
- **10.** Position the plastic ball in the second stage of the launcher.



STUDENT

- **11.** By pressing START on the display it will show the "standby" signal (*) flashing. By pressing START again, the signal (*) is in the operating mode, waiting for the ball to pass through the sensor.
- **12.** After the launch the screen shows the time interval value of the ball passing through the two sensors.
- **13.** Write down the time measured in the table 1. Repeat data collection at least 5 times.
- **14.** By pressing START the time measurement of the screen is sent to memory. For a new measure go back to item 11. Complete the table 1.

Table 1							
N	Δx (m)	Δt (s)	v (m/s)				
1	0,05						
2	0,05						
3	0,05						
4	0,05						
5	0,05						

15. Use the equation below and calculate the ball velocity through the sensors:

$$v = \frac{\Delta x}{\Delta t}$$

- **16.** Determine the average value of the launch velocity.
- **17.** Compare the experimental values of the launch velocities found by processes (a) and (b). Was there any significant difference between the values found?

Method	By the reach measure and launch height (v_{01})	By the time of passage between two sensors (v ₀₂)
Velocity (m/s)		, , , , , , , , , , , , , , , , , , ,

18. Calculate the average launch velocity of the two procedures.

19. Calculate the percentage error: $e^{0}\!\!/_{0}=rac{|v_{01}-v_{02}|}{rac{v_{01}+v_{02}}{2}}$


20. Justify possible inconsistencies between results.

Part III - Determination of horizontal velocity of launch (v_o) using the measure of the reach (A) and the time of the projectile in the air



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the cannon as shown in the figure (with a sensor and the flight time sensor). Cover the flight sensor with a carbon paper to mark the point where the ball touches the lowest point of the path.



- **2.** Connect the sensor and platform to the timer.
- **3.** Adjust the sensor position as close as possible to the launcher's mouth.
- **4.** Turn on the timer and check if it has identified the sensors and the platform. Test the sensors by interrupting the infrared passage. An intermittent exclamation mark (!) Should appear in the upper right corner of the screen indicating that the sensor 1 are identified. Tap the sensor briefly, the display shows signal 2 (!) indicating that the timer has identified the flight sensor.
 - When you turn on the timer, it displays the STANDBY screen and should show:
 - STB
 - Vb: 8,5V Output Voltage.
 - Sns :2- Number of sensors identified.
 - Bobbin: off indicating the coil is off.
- **5.** Select function 1 by clicking FUNC. The display will show the function 1 F1. The display should read "[F1]" "Btn nS:2". In this function the timer will be using two sensors, and the timer will present on time measurement. (See information in the timer manual).





- **6.** Schedule the timer for data acquisition by starting the time counting on the first sensor and interrupting the second one.
- 7. Press SETUP, press for a few seconds until the display shows CFG, use the key (↓) and adjust from Btn to Sns. Click START to end programming. On the nS: 1 screen, it indicates that the timer will only take one measure of time.



- **8.** Position the plastic ball in the second stage of the cannon.
- **9.** By pressing START on the display it will show the "standby" signal (*) flashing. By pressing START again, the signal (*) is in the operating mode, waiting for the ball to pass through the sensor.
- **10.** After the ball touches the flight sensor, the time counting stops.
- **11.** Write down in the table the measured time. Repeat the procedure at least 5 times and calculate the average reach distance (A). To measure the reach, observe the following precautions:
 - Place a plumb bob from the center of the sensor to the floor. Mark this reference.
 - The reach shall be measured from the reference point on the floor to the point where the ball touches the platform.
- **12.** By pressing START the time measurement of the screen is sent to memory. For a new measure go back to item 11. Complete the table 1.

 Table 1

 N
 A (m)
 Δt (s)
 v (m/s)

 1
 ...

 2
 ...

 3
 ...

 4
 ...

 5
 ...

 Average value



STUDENT	

13. Use the follow equation and calculate the ball velocity through the sensors:

$$\boldsymbol{v} = \frac{\boldsymbol{A}}{\Delta \boldsymbol{t}}$$

- 14. Determine the average value of the launch velocity.
- **15.** Compare the experimental values of the launch velocities found by processes (a) and (b). Was there any significant difference between the values found?

Method	By the reach measure and launch height (v ₀₁)	By the time of passage between two sensors (Vo2)	By the residence time in the air (v ₀ 3)
Velocity (m/s)			

16. Calculate the average value of the launch velocity found by the two procedures.

EXPERIMENT 02 - PROJECTILE OBLIQUE LAUNCH

Objectives:

- Recognizing the physical quantities involved in a projectile oblique launch.
- Check the relationship among the physical quantities present in an oblique launch.

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING
XX	62001226	01	Un.	DIGITAL TIMER AZB-30 USB (*)
XX	62001201	02	Un.	PHOTOELECTRIC SENSOR PGS-D10 (*)
XX	04002037	01	Un.	FLIGHT TIME SENSOR TFS-D10


(*) It does not accompany the product. It is sold separately.

Part I - Measure of Reach (A) and determining of launch velocity (v_0) using the measurement of the passage time between two sensors.



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the equipment as shown (the space between sensors is 5 cm):





S	T	U	D	E	N	Т						

- 2. Align sensors correctly and position the first of them as close as possible to the launcher's exit.
- **3.** Connect the sensors to the timer.
- **4.** Measure the distance between the two sensors centers.
- **5.** Adjust the cannon launch position to an angle of 30°.
- **6.** Measure the launch height (y_0) (from the bottom of the ball to the floor).
- 7. Launch the ball once to put a paper sheet and a carbon paper over it.
- **8.** Select function 1 by clicking FUNC. The display will show the function 1 F1. The display should read "[F1]" "Btn nS:2". In this function the timer will be using two sensors, and the timer will present two time measurements. (See information in the timer manual).



- **9.** Schedule the timer for data acquisition by starting the time counting on the first sensor and interrupting the second one.
- **10.** Press SETUP, press for a few seconds until the display shows CFG, use the key (↓) and adjust from Btn to Sns. Click START to end programming. On the nS: 1 screen, it indicates that the timer will only take one measure of time.



- **11.** Position the plastic ball in the second stage of the launcher.
- **12.** By pressing START on the display it will show the "standby" signal (*) flashing. By pressing START again, the signal (*) is in the operating mode, waiting for the ball to pass through the sensor.



STUDENT	

- **13.** Fire the launcher. After the launch the screen shows the time interval value of the ball passing through the two sensors.
- **14.** Measure the horizontal range of the ball (A) and the time interval value and note this data in the table. Repeat data collection at least 5 times.
- **15.** By pressing START the time measurement of the screen is sent to memory. For a new measure go back to item 12. Complete the table 1.

	Table 1					
N	A(m)	t _p (s)	v _o (m/s)			
1						
2						
3						
4						
5						
Average value						

16. Use the values in the table and calculate the average value of the measured reach - Ameasured

$$A_{med} = \frac{\sum A_i}{N}$$

17. Use the equation below and table values to calculate the ball velocity through the sensors:

$$\boldsymbol{v} = \frac{\Delta \boldsymbol{x}}{\Delta \boldsymbol{t}}$$

18. Determine the average value of the launch velocity (v_o).

19. Use the values obtained for the initial velocity (v_0) , height (y_0) and launch angle (θ) and determine the predicted reach value using the equation:

$$A = v_0 \cdot cos\theta \cdot \frac{v_0 \cdot sen\theta + \sqrt{v_0^2 \cdot sen^2\theta + 2gy_0}}{g}$$

STUDENT	

20. Compare the predicted reach value with the measured value.

$$e\% = \frac{\left|A_{measured} - A_{predicted}\right|}{\frac{A_{measured} + A_{predicted}}{2}}$$




Part II – Determining launch velocity (v_0) using the measure of reach (A) and the projectile residence time interval in the air.



Do not look directly at the cannon exit, as it may be charged!

1. Assembly as shown in the figure (with a sensor and the flight time sensor).



- **2.** Connect the sensors and the platform to the timer.
- **3.** Adjust the sensor position as close as possible to the launcher's mouth.
- **4.** Move the assembly to a launch at an angle of 60°.
- **5.** Turn on the timer and check if it has identified the sensors and the platform. Test the sensors by interrupting the infrared passage. An intermittent exclamation mark (!) should appear in the upper right corner of the screen indicating that the sensor 1 are identified. Tap the sensor briefly; the display shows signal 2 (!) indicating that the timer has identified the flight sensor.
 - When you turn on the timer, it displays the STANDBY screen and should show:
 - STB
 - Vb: 8,5V Output Voltage.
 - Sns :2- Number of sensors identified.
 - Bobbin: off indicating the coil is off.
- **6.** Select function 1 by clicking FUNC. The display will show the function 1 F1. The display should read "[F1]" "Btn nS:2". In this function the timer will be using two sensors, and the chronometer will present two time measurements. (See information in the timer manual).





- **7.** Schedule the timer for data acquisition by starting the time counting on the first sensor and interrupting the second one.
- **8.** Press SETUP, press for a few seconds until the display shows CFG, use the key (↓) and adjust from Btn to Sns. Click START to end programming. On the nS: 1 screen, it indicates that the timer will only take one measure of time.



- **9.** Position the plastic ball in the second stage of the launcher and perform some launches to identify the probable position of the platform.
- **10.** By pressing START on the display it will show the "standby" signal (*) flashing. By pressing START again, the signal (*) is in the operating mode, waiting for the ball to pass through the sensor.
- **11.** After the ball touches the flight sensor, the time counting stops.
- **12.** Note the time interval value provided by the timer in the table.
- **13.** Measure the reach (A) and note it in the table. To measure the reach, observe the following precautions:
 - Place a plumb bob from the center of the sensor to the floor. Mark this reference.
 - The reach shall be measured from the reference point on the floor to the point where the ball touches the platform.
- 14. Repeat data collection at least 5 times.



STUDENT	•
	STUDENT

15. By pressing START the time measurement of the screen is sent to memory. For a new measure go back to item 10.

N	A (m)	Δt (s)	v _{ox} (m/s)	v₀(m/s)
1				
2				
3				
4				
5				
	P	_		

16. Use the follow equation and calculate the ball horizontal velocity (v_{ox}) :

$$v_{ox} = \frac{A}{\Delta t}$$

17. Determine the average value of the launch velocity (v_0) for each measure carried out.

$$v_o = \frac{v_{ox}}{\cos 60^o}$$

- **18.** Determine the average value of the launch velocity (v_{om}) .
- **19.** Compare the experimental values of the launch velocities found by processes (a) and (b). Was there any significant difference between the values found?

	By the time of passage between two sensors (v ₀₂)	
Velocity (m/s)		

20. Compare the results.

EXPERIMENT 03 - RELATIONSHIP BETWEEN REACH AND LAUNCH ANGLE

Objective:

- Checking the dependence between a projectile's reach and launch angle.

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING


(*) It does not accompany the product. It is sold separately.

Part I – Launch without gap



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the launcher as shown.



- **2.** Provide a platform (for example, a box) that allows the projectile to reach it at the same horizontal point as the launch point.
- 3. Adjust the launch angle to 10° and insert the plastic ball into the second stage.
- **4.** Make some previous launches to locate the point where the ball will touch the surface of the platform.



STUDENT	

- **5.** Fix a paper sheet on the platform and a carbon paper over
- **6.** Make three launches and measure the value of launch position A (m) to the point where the ball touches the platform.
- **7.** Repeat the launch procedures for the angles suggested in table 1.

Table 1


Angle	Reach 1	Reach 2	Reach 3	Reach measured	Reach calculated
10					
20					
30					
40					
45					
50					
60					
70					
80					

- **8.** Use the experimental values from table 1 and calculate the average reach distance for each launch angle.
- **9.** Use launch velocity equal to 4.12 in expression:

$$A_{calc} = \frac{v_0^2 \times sen2\theta}{g}$$

And obtain the calculated reach value for each angle.

10. Draw the reaches graph (A_{measured}) according to the launch angle.



- **11.** According to the table and the graph, what is the launch angle value that provides the longest reach?
- **12.** What can be concluded about the reach value for the complementary angles, that is, 10° and 80°, 20° and 70°, 30° and 60°, 40° and 50°?

STUDENT

13. Are the experimental results consistent with the theory?




Part II - Launch with gap



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the launcher as shown so that the projectile reaches the floor.



- **2.** Use a plumb bob to determine the initial launch position on the floor.
- **3.** Adjust the launch angle to 10° and insert the plastic ball into the second stage.
- **4.** Fire at least three launches and measure the value of launch position A (m) to the point where the ball touches the floor. Fix a sheet of white paper on the platform and onto it a sheet of carbon paper.
- **5.** Fire at least three launches and measure the value of launch position A (m) to the point where the ball touches the floor.
- **6.** Repeat the launch procedures for the angles suggested in table 2.

	T	Table 2		
Angle	Reach 1 (m)	Reach 2 (m)	Reach 3 (m)	$A_{average}$
10°				
20°				
30°				
40°				
50°				
60°				
70°				
80°				



- **7.** Use the experimental values from table 2 and calculate the average reach value (A_{medido}) for each launch angle.
- 8. Use the equation:

$$A_{calculado} = v_0 \cdot cos\theta \cdot \frac{v_0 \cdot sen\theta + \sqrt{v_0^2 \cdot sen^2\theta + 2gy_0}}{g}$$

to obtain the reach calculated value


Table 3

	Tubic 3	
Angle	Acalculated (m)	A _{measured} (m)
10°		
20º		
30⁰		
40°		
50∘		
60⁰		
70°		
80⁰		

9. For each of the launch angles, are the reach values (measured and calculated) compatible?

	STUDENT	

10. Draw the reaches graph (A_{measured}) according to the launch angle.



- **11.** According to the table and the graph, what is the launch angle value that provides the longest reach?
- 12. Are the measured and calculated values in agreement with the maximum reach value?

EXPERIMENT 04 - ANALYSIS OF A TRAJECTORY OF A PROJECTILE

Objective:


- Obtaining the relationship between the vertical and horizontal position of the projectile during its motion.


Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62001074	01	Un.	STEEL BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING

Part I - Horizontal Launch



Do not look directly at the cannon exit, as it may be charged!





- 1. Assemble the launcher as shown.
- **2.** Adjust the cannon for horizontal launches. Use a plumb bob to determine on the floor the horizontal position of the projectile $(x_0 = 0)$.
- **3.** Assemble near a wall or set up a vertical bulkhead so that the ball, once launched, collides with it and determines the vertical and horizontal positions of the projectile in the collision with this target, as shown in the figure.




STUDENT	

- **4.** Position the bulkhead at x = 0.200 m from the horizontal launch position marked on the floor.
- **5.** Fix a paper sheet on the target place on the wall and with a carbon paper over it.
- **6.** Fire the ball and measure the value of the vertical position in which the projectile hits the target.
- **7.** Reposition the target to the horizontal positions suggested in the table and repeat the launch procedures until you complete it. Complete the table 1.

	Table 1	
X (m)	Y (m)	X^{2} (m ²)
0,000		
0,200		
0,400		
0,600		
0,800		
1,000		
1,200		
1,400		
1,600		
1,800		




8. Use the values collected in table 2 and draw the graph "vertical position (y)" versus "horizontal position (x)".



9. What is the aspect of the curve obtained?

10. Do the appropriate change of variables and linearize the graph.




- **11.** Obtain the equation that relates the quantities y x
- **12.** Combine the equations of the two movements performed by the ball and obtain analytically the equation that relates the two positions (y and x) of the ball in its trajectory.

Part II - Oblique launch



Do not look directly at the cannon exit, as it may be charged!

- **1.** Assemble the launcher as shown.
- 2. Adjust the launch angle by 60° and measure the height (yo) of the launch point in relation to the tabletop.
- **3.** Assemble on the table a vertical bulkhead so that the ball, once launched, collides with it and determines the vertical and horizontal positions of the projectile in the collision with this target, as shown in the figure.




- **4.** Use a plumb bob to obtain, on the tabletop, the initial horizontal position $(x_0=0)$.
- **5.** Use the steel ball and insert it into the second stage of the cannon.
- **6.** Position the bulkhead at 0,050 m from the horizontal launch position $(x_0 = 0)$ marked on the table top.
- **7.** Affix a paper sheet on the target place with a carbon paper over it.



8. Reposition the target to the horizontal positions suggested in the table and repeat the launch procedures until you complete it. Complete the table 1.

N	Table 1 x (m)	y (m)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		

9. Draw the graph "vertical position y" versus "horizontal position x".



10. Use the software feature (Excel) and obtain the equation corresponding to the curve.

STUDENT	

11.	Combine the equations of the motions	(vertical and horizontal)	and obtain the equation y =
	f(x).		

12. Is the dependence between y and x theoretically found consistent with the equation graphically found?

EXPERIMENT 05 - ANALYSIS OF A TRAJECTORY OF A PROJECTILE II

Objective:

- Verifying the validity of the equation that gives the projectile position at a given moment in an oblique launch.


Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING
XX	62001226	01	Un.	DIGITAL TIMER AZB-30 USB (*)
XX	62001201	01	Un.	PHOTOELECTRIC SENSOR PGS-D10 (*)
XX	04002037	01	Un.	FLIGHT TIME SENSOR TFS-D10
XX	29003002	01	Un.	BIG TRIPOD (*)
XX	30002014	01	Un.	PAIR OF ADAPTABLE STEMS 12,7MM X 400MM (*)
XX	31003001	01	Un.	METALLIC CLIP AZB-027

(*) It does not accompany the product. It is sold separately.



Do not look directly at the cannon exit, as it may be charged!

- **1.** Assemble the launcher at an angle of 40 ° to the horizontal on the side of the tabletop, so the projectile has space to move and fall on the floor, as shown in the figure.
- **2.** Assemble the flight time sensor as shown.



3. Using the plumb bob, affix a paper sheet on the floor aligning with the cannon mouth.



STUDENT	


- **4.** Measure the launch height (yo) in relation to the floor. The measurement shall be made from the center of the ball indicated at the exit of the cannon, to the floor.
- **5.** Position the flight time sensor at x = 0.400m from the launcher's mouth. Cover the face of the sensor with a carbon paper.
- 6. Place the plastic ball in the projectile launcher and compress the spring to the first stage. Pull the trigger and measure the height y (vertical position) in which the ball touches the platform. Repeat this procedure three times and note in the table the average value of y.
- **7.** Reposition the flight time sensor to the positions suggested in the table and repeat the experimental procedures to obtain the values of y.

N	Horizontal position x (m)	Vertical position y (m)
1	0,000	
2	0,200	
3	0,400	
4	0,600	
5	0,800	
6	1,000	
7	1,200	
8	1,400	
9	1,600	
10	1,790	



8. Combine the oblique launch equations of a projectile and obtain the expression that gives the vertical position y according to the horizontal position x, that is: y = f(x).

9. With the experimental data of y and x draw the graph y versus x.



- **10.** What does the graph y = f(x) look like?
- **11.** Obtain the equation that represents the curve obtained in the graph.

CTUDENT	
STUDENT	

12. Does the equation obtained experimentally agree with the theoretical expression?




EXPERIMENT 06 - TRAJECTORY OF A PROJECTILE IN THE HORIZONTAL LAUNCH

Objective: Verifying if a horizontally launch projectile describes a parabolic trajectory.

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING
XX	29003002	04	Un.	BIG TRIPOD (*)
XX	30002014	06	Un.	PAIR OF ADAPTABLE STEMS 12,7MM X 400MM (*)
XX	02005006	04	Un.	LAB CLAMP WITH RING (*)
XX	31003001	04	Un.	METALLIC CLIP AZB-02 (*)

(*) It does not accompany the product. It is **s**old separately.

13. Assemble the launcher as shown.



1. Measure the launch height (yo) in relation to the floor. The measurement shall be made from the lower part of the ball indicated at the exit of the cannon, to the floor, according to the figure.



2. Affix a paper sheet and mark the vertical alignment of the plumb with the center of the ball, as shown.



- **3.** Place the plastic ball in the cannon and compress the spring to the first stage. Pull the trigger and observe where the ball touches the floor.
- **4.** At the dropping point of the ball, affix a paper sheet and on it a carbon paper, marking the projectile's reach distance A.
- **5.** Repeat the launch 5 times and measure the reached distance A by using the measuring tape.

Table 1				
N	Reach (m)			
1				
2				
3				
4				
5				
Average value				

6. Determine the average reach distance.

7. By combining equations:

$$\begin{cases} A = v_0 \cdot t \\ y = y_0 - \frac{g \cdot t^2}{2} \end{cases}$$

We obtain the equation that provides the launch velocity:

$$v_0 = A \cdot \sqrt{\frac{g}{2y_0}}$$

- **8.** Calculate the projectile horizontal launch velocity v_0 by using the reach average value.
- **9.** Draw on the sheet of paper placed on the floor a line that joins the point marked with the plumb bob and the point corresponding to the reach average value.
- **10.** Mark from the origin along the line drawn on the paper the following positions:

$$x_1 = 0.30 \text{ m}$$
; $x_2 = 0.60 \text{ m}$; $x_3 = 0.90 \text{ m}$; $x_4 = 1.20 \text{ m}$.

11. The vertical motion equation is:

$$y = y_o - \frac{g}{2}t^2; \ (v_{0y} = 0)$$

By combining with horizontal motion equation, $x = v_0 \cdot t$, results the equation that relates the projectile positions y and x, y=f(x):

$$y = y_0 - \frac{g}{2v_0^2}x^2$$

12. Calculate the value of the vertical y position of the projectile for x positions marked on the sheet of paper and suggested in the table.

Table 2

x (m) y₀ (m) v₀ (m/s) y (m)

0,300

0,600

0,900

1,200

13. Assemble four supports for the rings and position them in x positions suggested in table 2, as shown in the figure.



14. Prepare the projectile launch with the plastic ball in the first stage. Fire the launch and verify if the projectile has passed through the rings without touching them.

EXPERIMENT 07 - CONSERVATION OF ENERGY

Objective:

Verifying conservation of mechanical energy using the vertical launch.


Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
04	62001074	01	Un.	STEEL BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING
XX	62001226	01	Un.	DIGITAL TIMER AZB-30 USB (*)
XX	62001201	01	Un.	PHOTOELECTRIC SENSOR PGS-D10 (*)

(*) It does not accompany the product. It is sold separately.



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the projectile launcher as shown.



- 2. Adjust the launch angle to 90°.
- **3.** Measure the diameter (\emptyset) and the mass (m) of the steel ball.



- **4.** Assemble a tripod with a stem on the table and fix an object (a ruler, for example) that serves as a reference to determine the vertical reach of the projectile, as shown in the figure.
- **5.** Fix the sensor to measure the ball passage time, in the position shown in the figure.
- **6.** Connect the timer and adjust it to use the function 2, with the screens shown below:





To select the desired measurement type, press the SETUP key for 2s to enter the configuration mode. In the configuration mode, use the START / RESET key to navigate among the configuration parameters. To change a selected parameter use the SETUP / MEM \uparrow keys. After setting the function, press the FUNC key to save the selected parameters.

- **7.** Position the steel ball in the second stage of the launcher and fire the cannon 5 times to determine the vertical reach value of the ball, adjusting the referential at each launch.
- **8.** Note in table 1 the vertical reach and the time of passage values of the ball by the sensor.

Table 1					
N	y (m)	t(s)			
1					
2					
3					
4					
5					
y average					



CTUDENT
STUDENT

9.	Calculate the average value of the maximum height (y) reached by the ball.
10.	Calculate the average value of the passage time (t).
11.	Use the ball diameter value and the average passage time and calculate the launch velocity modulus $\nu_{0}. \label{eq:constraint}$
12.	Calculate the kinetic energy value in the initial position (at the time of launch).
13.	Calculate the value of the gravitational potential energy at the highest point of the path.
14.	Compare the values of mechanical energy in both procedures and justify the discrepance found:
15.	Does the experiment confirm the principle of conservation of mechanical energy?

EXPERIMENT 08 - CONSERVATION OF THE LINEAR MOMENT

Objective:

Verifying the conservation of linear momentum in elastic and inelastic collision.


Item	Code	Quant.	Unit	Description			
01	62002176	01	Un.	CLAMP "C"			
02	62005751	01	Un.	CANNON LAUNCHER HOLDER			
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)			
04	62001074	02	Un.	STEEL BALL Ø25MM			
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION			
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)			
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)			
09	62005317	01	Un.	CANNON'S FASTENER			
10	62001023	01	Un.	SENSORS/SPHERE SHOCK'S HOLDER			
11	62002056	01	Un.	BALL'S MAGNETIC FIXING (FOR COLLISION)			
12	53003001	01	Un.	FASTENER FOR BALL'S MAGNETIC FIXING			
13	03003011	01	Un.	TAPE MEASURE 05M			
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING			

Part I - Elastic Oblique Collision




Do not look directly at the cannon exit, as it may be charged!

1. Assemble the projectile launcher as shown.






- 2. Position the launcher for horizontal launch
- **3.** Measure the mass of each steel ball.
- **4.** Attach the magnetic fixing as shown previously.
- **5.** Fire the cannon once with the steel ball 1 in the second stage and adjust the position of the launcher to ensure that the ball falls on the sheet of paper covering the table.
- 6. Place the second steel ball (2) on the magnetic fixing and adjust its position very well for an oblique collision at an angle Θ around 40 ° to 50 °. It is advisable that the collision occurs at least 3.0 cm from the launcher's mouth.



- 7. Fire launches that produce a collision between steel balls 1 and 2 and make adjustments to the launcher and the support of ball 2 positions that ensure that in these procedures (with and without collision) the balls fall on the table.
- **8.** Use a plumb bob and determine on the sheet of paper placed on the table the point O where the collision occurs. This point will be considered as the origin of horizontal displacement.



- **9.** Measure the launch height h_0 in relation to the table surface.
- **10.** Remove the ball 2 from the magnetic support.
- **11.** Cover the sheet of paper with carbon paper in the estimated position.
- **12.** Fire five launches with ball 1 in the first stage and determine the point where the ball touches the table surface.



STUDENT	

13. Measure the reach value in each launch, note them in table 1 and determine the average value A_{\circ} .


Tale 1				
N	Height ho (m)	Reach A _o (m)		
1				
2				
3				
4				
5				
Average val	ue			

- **14.** Place the ball 1 in the first stage of the launcher. Position the ball 2 on the magnetic fixing and fire the launcher. Observe where the balls touch the table surface after the collision. Place in these positions sheets of carbon paper.
- **15.** Prepare a new launch for the two balls collision.
- **16.** Fire the launcher and measure the two balls reaches A_1 and A_2 after the collision.
- **17.** Repeat the launch process five times and determine the average position of each ball reach after the collision. Complete table 2

Table 2					
event	A_1 (m)	$A_2(m)$			
1					
2					
3					
4					
5					
Average					
value					

- **18.** Draw on the sheet of paper from the point O, the lines corresponding to the average reaches A₀ of the ball 1 before the collision, A1 of the ball 1 and A2 of the ball 2, after the collision.
- 19. Measure the angle Θ between the vectors corresponding to A1 and A2.






20. Use the average value of table 1 and the appropriate equations to calculate the initial momentum modulus p_0 of the system:

$$v_o = A \cdot \sqrt{\frac{g}{2h_o}} e p_o = m_{esf} \cdot v_o$$

- **21.** Use the data in table 2 and calculate the values of linear momentum modules of each ball after the collision.
- **22.** Scale the vectors \vec{p}_o , \vec{p}_1 and \vec{p}_2 .



- **23.** Apply the cosine law and obtain the vector module \vec{p} , resulting from the vectors \vec{p}_1 and \vec{p}_2 .
- **24.** Compare the values of the linear momentum modules just before and immediately after the collision. Has the linear momentum of the system been conserved?

Table 3

Initial linear	Final linear	
momentum p₀	momentum	Percent error
(kg.m/s)	(kg.m/s)	

- **25.** Take the direction x of reference as the direction of the ball's reach before the collision and trace the pair of axes (x, y). Measure the angles that each vector forms with the x-axis.
- **26.** Obtain the orthogonal components of the linear moments after the collision:

Table 4

		Angle	with	Component x	Component
	Module	x-axis		X	У
				p_{ox}	p_{oy}
p_o					
p_1					
p ₂					
р					

27. Can we consider that in the direction x the linear momentum has been conserved? And in the direction y?

28. Was there conservation of kinetic energy? Can collision be considered as an elastic collision? Justify.

Table 5

	Final ki	Final kinetic energy (J)			
Initial kir	netic				
energy (J)	Ball 1	Ball 2	Sum	e%	

Part II - Inelastic Oblique Collision



Do not look directly at the cannon exit, as it may be charged!

- 1. Use the same assembly of the first part with steel balls.
- 2. Consider the data from the table obtained with the ball 1 launch without collision.

	Table 1	
N	Launch height ho (m)	Reach A _o (m)
1		
2		
3		
4		
5		
A۱	erage value	


- **3.** Adopt as axis x the launch direction without collision.
- **4.** Wrap the ball 2 with one lap of sticky tape and repeat the collision launch procedures. In this case it is advisable to Fire only one launch, since the replacement of the ball 2 would hardly occur in the same positioning conditions
- **5.** Lay on the table a sheet of paper where the balls fall after the collision.
- **6.** Repeat the procedures for marking the point of origin of the horizontal displacements of the balls.
- **7.** Fire some previous launches and observe the points where the balls touch the table. Put sheets of carbon paper at these points.
- 8. Fire the cannon and score the points obtained from the two balls reaches.

	Reach (m)	Angle with axis x
Ball 1		
Ball 2		

9. Calculate the value of the linear momentum module of the system immediately before the collision.



10. Draw on the sheet of paper from the point O, the lines corresponding to the average reaches A0 of the ball 1 before the collision, A1 of the ball 1 and A2 of the ball 2, after the collision.



- **11.** Measure the angle that each vector forms with the vector direction obtained without the collision.
- **12.** Measure the angle Θ between the vectors corresponding to A1 and A2.
- **13.** Calculate the value of the velocity and momentum modules of each ball after the collision. The launch height h_0 is the same as the first part. ($h_0 = 0,189$ m).

$$v = A \cdot \sqrt{\frac{g}{2h_o}} ep_o = m_{esf} \cdot v_o$$

- 14. Use the cosine law and calculate the value of the linear momentum module after the collision.
- 15. Was there conservation of the linear momentum of the system in the inelastic collision?

STUDENT	SIOPLIII
---------	----------

16. Calculate the initial and final kinetic energy of the system. Was there energy conservation in the collision?



Part III - Frontal Elastic Collision



Do not look directly at the cannon exit, as it may be charged!

- **1.** Use the same assembly as in the first part.
- 2. Consider the data of table 1 obtained with the ball 1 launch without collision.

N	Table 1 Launch height h ₀ (m)	Reach A _o (m)
1		
2		
3		
4		
5		
Averag	e value	

- **3.** Adopt as x-axis the launch direction without collision.
- **4.** Adjust the position of the ball 2 support so the centers of the two balls are on the same axis providing a frontal collision.
- **5.** Lay on the table a sheet of paper to accommodate the balls after the collision.
- **6.** Fire the cannon a couple of times and adjust the position of ball 2 until a really frontal collision (without lateral deviation of the balls after the collision).
- **7.** Cover the sheet with carbon paper in the regions where the balls touch the table.
- 8. Trigger the launcher and mark the reaches' points A1 and A2 of the two balls.

9. Draw the lines corresponding to the scales, measure and note the value of the reaches modules A_1 and A_2 of each ball.



	STUDENT	
--	----------------	--

10. Measure the angles Θ_1 and Θ_2 which A_1 and A_2 form with the x-axis and the angle Θ between A_1 and A_2 .

	Table 2		
	Reach	Angle	with
	(m)	x-axis	
Ball 1			
Ball 2			

11. Calculate the value of the velocity and momentum modules of each ball after the collision. The launch height h_0 is the same as the first part. ($h_0 = 0,189$ m).

$$v_o = A \cdot \sqrt{\frac{g}{2h_o}} e p_o = m_{esf} \cdot v_o$$

- **12.** Obtain the linear momentum module after the collision.
- 13. Was there conservation of the linear momentum of the system in the inelastic collision?

14. Calculate the initial and final kinetic energy of the system. Was there energy conservation in the collision?

Part IV - Frontal Inelastic Collision



Do not look directly at the cannon exit, as it may be charged!

- **1.** Use the same assembly as in the first part.
- 2. Consider the data of table 1 obtained with the ball 1 launch without collision.

	Table 1			
N	Launch height h _o (m)	Reach A _o (m)		
1				
2				
3				
4				
5				
Avera	Average value			

- **3.** Adopt as x-axis the launch direction without collision.
- **4.** Wrap ball 2 with one lap of sticky tape and repeat collision launch procedures. In this case it is advisable to fire only one time, since the replacement of the ball 2 would hardly occur in the same positioning conditions.
- **5.** Adjust the position of the ball 2 support so the centers of the two balls are on the same axis providing a frontal collision.
- **6.** Lay on the table a sheet of paper to accommodate the balls after the collision.
- **7.** Fire a couple of times and adjust the position of ball 2 until a really frontal collision (without lateral deviation of the balls after the collision).
- **8.** Cover the sheet with carbon paper in the regions where the balls touch the table.
- **9.** Fire the launcher and mark the reaches' points A1 and A2 of the two balls.

- **10.** Draw the lines corresponding to the scales, measure and note the value of the reaches modules A1 and A2 of each ball.
- **11.** Measure the angles Θ_1 and Θ_2 which A_1 and A_2 form with the x-axis and the angle Θ between A_1 and A_2 .



Table 2

	Tubic 2	
	Reach	Angle with x-
	(m)	axis
Ball 1		
Ball 2		

12. Calculate the value of the velocity and momentum modules of each ball after the collision. The launch height h_0 is the same as the first part. ($h_0 = 0,189$ m).

$$v_o = A \cdot \sqrt{\frac{g}{2h_o}} e p_o = m_{esf} \cdot v_o$$

- **13.** Obtain the linear momentum module after the collision.
- 14. Was there conservation of the linear momentum of the system in the inelastic collision?

15. Calculate the initial and final kinetic energy of the system. Was there energy conservation in the collision?

EXPERIMENT 09 – LAUNCH ANGLE THAT MAXIMIZES THE HEIGHT OF A PROJECTILE

Objective:


- Finding the launch angle that maximizes the height reached by projectile launched from a fixed distance from a bulkhead.

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
05	62005120	01	Un.	PLASTIC BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the launcher as shown.



- 2. Position the table facing a wall so the launcher's mouth is about 1.60 m away. Use a plumb bob to determine the initial launch position on the floor and measure the distance (x) from the wall to the launcher's mouth and the vertical launch position (y_0) .
- **3.** Fire the cannon and determine where to place a 2 sheet of paper, displayed vertically.



STUDENT	

- **4.** Adjust the launch angle to 10° and insert the plastic ball into the second stage.
- **5.** Fire launches to locate the point where the ball touches the wall. Attach at this point a sheet of carbon paper.
- **6.** Fire the launch and measure the height (from the floor) reached by the projectile.
- **7.** Launch at the angles suggested in table 1 and note the respective heights reached.

					Ta	ble 1						
Angle θ	10	15	20	25	30	35	40	45	50	55	60	65
(degree))											
Height												
y (m)												

- **8.** Note in table 2 the launch velocity value obtained in experiment 1.
- **9.** Use the initial velocity obtained in experiment 1.
- **10.** By using the initial velocity of step 1 and the distance from the wall to the launcher, calculate the angle that provides the maximum height.



STUDENT

11. Draw the graph of the vertical position versus the launch angle of table 1.




Table 2

	Valor
Angle for maximum height-measured	
Maximum height	
Horizontal wall distance	
Launch height	
Initial velocity calculated	
Angle for maximum height-calculated	
Percentage difference among angles	

- **12.** Does the value found analytically agree with the values obtained experimentally? Justify.
- **13.** For the angle that gives the maximum height when the ball hits the wall, has it reached the peak of the trajectory?

STUDENT

14.	How far from the wall will the height be maximized for a launch at a 45° angle? What would
	be the maximum height in this case?

15. Launch to x equal to the value found for the 45° angle and measure the y value (x; 45°). Compare with calculated value.

ST	U	DE	EN	IT						

EXPERIMENT 10 - BALLISTIC PENDULUM - APPROXIMATE METHOD

Objectives:

- Using the conservation of linear momentum and conservation of mechanical energy in a ballistic pendulum to find the launch velocity of a projectile.
- Comparing the launch velocity found with the velocity obtained in a horizontal launch by using different processes.

PROJECTILE LAUNCHER

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
04	62001074	01	Un.	STEEL BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING

Ballistic Pendulum Accessories (sold separately)

Item	Code	Quant.	Unit	Description				
01	62005611	02	Un.	50G CYLINDER MASS				
03	53003001	01	Un.	PENDULUM'S TOWER THUMB SCREW'S NUT				
04	53001009	01	Un.	PENDULUM'S THUMB SCREW				
04	53004002	01	Un.	METAL HANDLE M3X10				
05	62002055	01	Un.	BALLISTIC PENDULUM				
06	62002182	01	Un.	PENDULUM'S TOWER				


Part I - Obtaining the horizontal launch velocity by using the reach measure.



Do not look directly at the cannon exit, as it may be charged!

1. Assemble the ballistic pendulum as shown.





- 2. Measure the launch height (yo) in relation to the floor. The measurement shall be made from the lower part of the ball indicated at the exit of the cannon, to the floor, according to the figure.
- **3.** With a plumb bob, mark on a sheet of A4 paper pasted with adhesive tape on the floor, the post position (origin of horizontal displacement). The plumb bob must match the vertical passing through the center of the ball.



- **4.** Place the steel ball in the projectile launcher and compress the spring to the third stage. Pull the trigger and observe where the ball touches the floor.
- **5.** At the dropping point of the ball, affix a paper and on it a sheet of carbon paper for marking the projectile reach distance A.



	STUDENT	

6. Repeat the launch 5 times and measure the reach distance A.

Table 1

N Reach (m)

1

2

3

4

5

Average value


7. Use the horizontal launch equations and calculate the launch velocity:

Part II - Obtaining launch velocity by using the conservation of linear momentum in a ballistic pendulum



Do not look directly at the cannon exit, as it may be charged!

- 1. Use the 1st part launcher assembly at an end of the table.
- 2. Measure the mass m of the metal ball and the mass M of the set (ballistic pendulum + ball). m=0.064 kg e M = 0.294 kg
- **3.** Insert the ball into the pendulum receiver and use a line to find the center of mass position of the set. Slide the line holding the pendulum until it remains in horizontal balance.
- **4.** Measure the distance R from the center of mass position to the axis of rotation.



- **5.** Attach the pendulum to the launcher support and make the necessary adjustments.
- **6.** For the correct positioning of the pendulum the following steps are suggested:
 - Keep the pendulum free in vertical.
 - Carefully approach the cannon until it softly engages the mouth of the pendulum trimmer.
 - Tighten the cylinder locking screws.


With these procedures the ball, when launched, will always be correctly picked up by the pendulum receiver.

7. Place the steel ball in the projectile launcher and compress the spring to the second stage.



STUDENT		

- **8.** Adjust the pendulum position so the axis of trimmer fits correctly into the launcher's mouth.
- **9.** Adjusting the zero of the angle marker.
- 10. Release the trigger for launch and note the angle value reached by the pendulum.
- **11.** Write down the measurements of angles of the pendulum.
- **12.** Choose the most repeatable measure and adjust the indicator to about four degrees less than that angle.
- **13.** Repeat the launch three times and calculate the average value of the angle.



The height h reached by the center of mass of the pendulum is calculated by:

$$h = R(1 - \cos\theta)$$

Considering that the linear momentum of the system is conserved both in an elastic collision and in an inelastic collision, we have:

$$p_{before\ collision} = p_{after\ collision} \rightarrow mv_0 = M \cdot v \quad (1)$$

After the collision the ballistic pendulum moves and its center of mass reaches a height h: As mechanical energy is conserved in this movement, we have:

$$K_{after\ collision} = U_{max\ hieght} \rightarrow \frac{Mv^2}{2} = MgR(1 - cos\theta) \rightarrow v = \sqrt{2 \cdot g \cdot R(1 - cos\theta)}$$
 (2)

By combining (1) and (2), we obtain the value of the launch velocity before the pendulum collision (v_0) :



$$v_o = \frac{M}{m} \sqrt{2 \cdot g \cdot R(1 - cos\theta)}$$

14. Use the deduced equation to calculate the ball launch velocity value.

15. Compare the result found with the average velocity of the ball obtained in the horizontal launch in the first part.

16. Calculate the percentage difference between the two results found: $d\% = \frac{|A-B|}{\frac{A+B}{2}} \times 100\%$

17. Does the percentage difference obtained confirm the validity of the conservation principle of linear momentum? Justify.

18. What can be concluded about the energy conservation in the collision?

STUDENT	

EXPERIMENT 11 - BALLISTIC PENDULUM-EXACT METHOD

Objectives:

- Using the conservation of angular momentum and conservation of mechanical energy in a ballistic pendulum to find the launch velocity of a projectile.
- Using the conservation of linear momentum and conservation of mechanical energy in a ballistic pendulum to find the launch velocity of a projectile.

PROJECTILE LAUNCHER

Item	Code	Quant.	Unit	Description
01	62002176	01	Un.	CLAMP "C"
02	62005751	01	Un.	CANNON LAUNCHER HOLDER
03	62002015	01	Un.	CANNON (PROJECTILE LAUNCHER)
04	62001074	01	Un.	STEEL BALL Ø25MM
06	62005177	01	Un.	TUBE FOR CANNON COMPRESSION
07	48005003	02	Un.	BUTTERFLY NUT (CANNON'S FASTENER)
08	50001004	02	Un.	FLAT WASHER (CANNON'S FASTENER)
09	62005317	01	Un.	CANNON'S FASTENER
13	03003011	01	Un.	TAPE MEASURE 05M
14	62005274	01	Un.	PLUMB BOB WITH MAGNETIC FIXING


Ballistic Pendulum Accessories (sold separately)

Item	Code	Quant.	Unit	Description
01	62005611	02	Un.	50G CYLINDER MASS
03	53003001	01	Un.	PENDULUM'S TOWER THUMB SCREW'S NUT
04	53001009	01	Un.	PENDULUM'S THUMB SCREW
04	53004002	01	Un.	METAL HANDLE M3X10
05	62002055	01	Un.	BALLISTIC PENDULUM
06	62002182	01	Un.	PENDULUM'S TOWER




Do not look directly at the cannon exit, as it may be charged!

- **1.** Assemble the launcher as in the previous experiment.
- **2.** Measure the mass m of the metal ball and the mass M of the set (ballistic pendulum + ball). $m = 0.064 \ kg \ M = 0.294 \ kg$
- **3.** Measure the distance r from the center of the ball to the axis of rotation of the pendulum. r = 0.255 m





- **4.** Insert the ball into the pendulum receiver and use a line to find the center of mass position of the set. Slide the line holding the pendulum until it remains in horizontal balance.
- **5.** Measure the distance R from the position of the center of mass to the axis of rotation.




- **6.** Attach the pendulum to the launcher support and make the necessary adjustments.
- **7.** For the correct positioning of the pendulum the following steps are suggested:
 - Keep the pendulum free in vertical.
 - Carefully approach the cannon until it softly engages the mouth of the pendulum trimmer.
 - Tighten the cylinder locking screws.

With these procedures the ball, when launched, will always be correctly picked up by the pendulum receiver.

- **8.** Place the steel ball in the projectile launcher and compress the spring to the second stage.
- **9.** Adjust the pendulum position to align correctly into the launcher's mouth.
- **10.** Adjusting the zero of the angle marker.
- 11. Trigger the cannon and note the angle value reached by the pendulum.
- **12.** Fire some measurements of the reach angle of the pendulum.
- **13.** Choose the most repeated measure and adjust the indicator to about four degrees less than that angle.



14. Repeat the launch three times and calculate the average value of the angle.



- **15.** Assemble the pendulum without the launcher so it can swing freely.
- **16.** Use a timer and measure the time of 20 complete oscillations of small amplitude. Repeat the procedure for at least three times and note the time value (t = 20T) for the 20 oscillations in the table.

The height h reached by the center of mass of the pendulum is as shown in (c):

$$h = R(1 - \cos\theta)$$

The potential energy stored in the system at the (c)

$$U = M \cdot g \cdot h = MgR(1 - cos\theta)$$

The kinetic energy immediately after the collision equals the potential energy at the point where the pendulum reaches the largest angle, and therefore:

$$U = K_d = M \cdot g \cdot h = MgR(1 - cos\theta)(1)$$

The kinetic energy K_d and the angular momentum L_d of the system immediately after the collision of the ball with the pendulum are given by the equations:

$$K_d = \frac{I \cdot \omega^2}{2}$$
 and $L_d = I_{conj} \cdot \omega$

By combining the two equations we can obtain the relation between L_d and K_d :

$$L_d = \sqrt{2 \cdot I_{conj} \cdot K_d} \qquad (2)$$

The angular momentum of the system immediately before the collision (L_a) is restricted to the angular momentum of the ball, since at that moment the pendulum is at rest:

$$L_a = I_{ball} \cdot \omega = m \cdot r^2 \omega = m \cdot r \cdot v_0 \quad (3)$$

$$\text{Where:} \begin{cases} m-\text{mass of the ball} \\ r-\text{distance from the ball to the axis of rotation} \\ v_0-\text{velocity of the ball before the collision} \end{cases}$$

Considering that the angular momentum L of the system is conserved and by combining (2) and (3):

$$L_{before\ collision} = L_{after\ collision}$$

$$m \cdot r \cdot v_0 = \sqrt{2 \cdot I_{set} \cdot K_d} \qquad (4)$$

By combining (1) and (4), we obtain:
$$v_0 = \frac{1}{mr} \sqrt{2 \cdot I_{set} \cdot M \cdot g \cdot R(1 - cos\theta)}$$
 (5)

In the angular displacement of the pendulum the torque: $au = I_{set} \cdot \alpha = -R \cdot M \cdot g \cdot sen\theta$

For small angles this equation can be written:
$$\alpha = \frac{d^2\theta}{dt^2} = -\frac{MgR}{I_{set}} \cdot \theta$$

Which is similar to a simple harmonic motion equation $\alpha = -\frac{k}{m}x = -\omega^2 x$

What allows us to write:
$$\omega^2 = \frac{MgR}{I} \rightarrow I_{conj} = \frac{MgR}{\omega^2} \rightarrow I_{conj} = \frac{MgRT^2}{4\pi^2}$$
 (6)

The period T can be obtained by making the pendulum perform small oscillations.

17. Calculate the average value of the period T of the pendulum.

Table 2										
Event	1	2	3	Average value						
Time of 20 oscillations (s)										
Period T (s)										

- 18. Use the equation (6) and calculate the moment of inertia of the set Iset.
- **19.** Use the equation (5) and calculate the ball velocity immediately before launch.
- **20.** Calculate the percentage difference between the two results found: $d\% = \frac{|A-B|}{\frac{A+B}{2}} \times 100\%$

S'	STUDENT

21. Does the percentage difference obtained confirm the validity of the conservation principle of angular momentum? Justify.

22. What can be concluded about the energy conservation in the collision?