

3. Apply to the dynamometer a tractive force approximately equal to 0,10N.
4. Turn on the equipment and adjust the oscillator to a frequency close to 30 Hz. Move the rod that fixes the dynamometer to obtain the first vibration mode (1st harmonic).

Note:

If the frequency or rating is too high, the security system of the device will be triggered and the message

5. Once the first harmonic is obtained, adjust the frequency to 30 Hz. Wait a few seconds for the frequency to be practically stable and note the power value indicated by the dynamometer, the number of nodes, the number of wombs and the wavelength.
 -
 -
 -
 -
6. Move the dynamometer up or down while simultaneously adjusting the frequency to 30 Hz to obtain the highest amplitude standing wave.




Setting the standing wave in the fundamental mode.

7. Repeat the procedures to collect the values of the tractive force (F) and the wavelength (λ) for the 2nd, 3rd and 4th modes of vibration, always maintaining the frequency at 30 Hz.

- Table 1


Vibration Mode	number of nodes	number of wombs	F(N)	$\lambda(m)$	$\lambda^2(m^2)$
1st					
2nd					
3rd					
4th					

9. Draw the graph $F=f(\lambda)$ (tractive force according to wavelength).



- 10.** What does the curve look like?
- 11.** What is the likely relationship between F and λ ?

- 12.** Use the appropriate variable change to linearize the graph.



- 13.** Obtain with the aid of Excel the equation that relates the strength and the wavelength, $F=f(\lambda)$ and the relation of proportionality between the two quantities.

Code	Quant.	Unit	Description
63002049	01	UN.	STATIONARY WAVE GENERATOR BASE WITH ROTATIONAL SYSTEM, ELECTRONIC CONTROL AND FREQUENCY METER
63005071	01	UN.	EXPANDABLE ROD FOR DYNAMOMETER
63005200	01	UN.	EXPANDABLE ROD WITH PULLEY
62001247	01	UN.	DYNAMOMETER 1N
38039005	01	UN.	POWER SUPPLY DC 12V/2A
48005004	01	UN.	BUTTERFLY NUT ¼"
53004003	01	UN.	THUMB SCREW
50002003	01	UN.	WASHER ¼"
63005192	01	UN.	ROPE DENSITY 3


The image shows a laboratory apparatus for measuring the surface tension of liquids. It consists of a black base, a vertical rod with a red and black handle, a white spherical weight, and a digital display showing '300'.

1. Use the same assembly as in the first part.
2. Fix the three-thread rope at the appropriate location and measure the rope length.
3. Apply to the dynamometer a tractive force approximately equal to 0,10N.
4. Turn on the equipment and adjust the oscillator to a frequency close to 10 Hz. Move the stem that fix the dynamometer to obtain the first vibration mode (1st harmonic), with the greatest possible amplitude of the standing wave.
5. Wait a few seconds for the frequency to be practically stable and note in the table the force value indicated by the dynamometer and the frequency of the wave.
6. Increase the frequency by approximately 3 to 4 Hz. Move the dynamometer in order to obtain the maximum amplitude of the standing wave corresponding to 1st harmonic and note the strength and frequency values.

7. Repeat the previous procedure for new frequency values always using the first harmonic (up to a maximum value of approximately 30 Hz).


[illegible][illegible]

- 8.** Draw the graph $F=f(\lambda)$ (tractive force according to wavelength).



- 9.** What does the curve look like?
- 10.** What is the likely relationship between F and f ?
- 11.** Use the appropriate variable change to linearize the graph.

Frequency f (Hz)	Strength F (N)	$f^2(\text{Hz})^2$




- 12.** Obtain with the aid of Excel the equation that relates the strength and the wavelength, $F=f(\lambda)$ and the relation of proportionality between the two quantities.

Code	Quant.	Unit	Description
63002049	01	UN.	STATIONARY WAVE GENERATOR BASE WITH ROTATIONAL SYSTEM, ELECTRONIC CONTROL AND FREQUENCY METER
63005074	01	UN.	EXPANDABLE ROD FOR DYNAMOMETER
	01	UN.	ROPES
56001001	10	UN.	HOOK TYPE "S"
62001247	01	UN.	DYNAMOMETER 1N
38039005	01	UN.	POWER SUPPLY DC 12V/2A


-

9. Use the required variable change and linearize the graph.



10. According to the linearized graph, what should be the relation between the frequency and wavelength of the standing wave in a rope subjected to a constant tractive force?
11. What is the physical meaning of the angular coefficient of the linearized graph?

Code	Quant.	Unit	Description
63002049	01	UN.	STATIONARY WAVE GENERATOR BASE WITH ROTATIONAL SYSTEM, ELECTRONIC CONTROL AND FREQUENCY METER
63005071	01	UN.	EXPANDABLE ROD FOR DYNAMOMETER
63005200	01	UN.	EXPANDABLE ROD WITH PULLEY
62001247	01	UN.	DYNAMOMETER 1N
38039005	01	UN.	POWER SUPPLY DC 12V/2A
48005004	01	UN.	BUTTERFLY NUT ¼"
53004003	01	UN.	THUMB SCREW
50002003	01	UN.	WASHER ¼"
63005194	01	UN.	COMBINED ROPE



1. Use the same assembly as in the third part.
2. Attach to the equipment the rope consisted of two parts with different densities.
3. Turn on the oscillator and adjust the frequency to approximately 25 Hz.
4. Adjust the dynamometer by moving it up or down until you find the first vibration mode, in the thinnest part of the rope.
5. Note the wave frequency value presented by the oscillator.

[illegible]

6. Consider the thinnest part of the rope with wavelength λ_1 and the thickest part with wavelength λ_2 .
7. Obtain the wavelengths of each wave propagating in the two means.
8. The wave propagates in the rope from the mean 1 (less dense) to mean 2 (denser).
9. When changing mean, did the wave frequency change? Justify.
10. What should be the relationship between the wave velocities in the two means?